

Bottles, square, wide neck, PP

NALGENE®

Square shape saves space. With leakproof polypropylene closure. Autoclavable.

Cat. No.	Code	Capacity, ml	Pk Qty
BO/36232	2110-0002	60	12
BO/36234	2110-0006	175	12
BO/36236	2110-0008	250	12
BO/36238	2110-0016	500	12
BO/36240	2110-0032	1000	6



BO/3623

Bottles, heavy duty, wide neck, PP, Autoclavable

NALGENE®

Large, thick-walled bottles designed for heavy-duty applications. Robust polypropylene closure with TPE (thermoplastic elastomer) gasket guarantees leakproof service under vacuum. Also useful as waste aspirator bottle.

Cat. No.	Code	Capacity, litres	Pk Qty
BO/36242	2126-1000	1	6
BO/36244	2126-2000	2	2
BO/36246	2126-4000	4	1



BO/3624

Bottles, carboy, wide neck, PP, Leakproof

NALGENE®

Popular round shape carboy with handles and white polypropylene leakproof screw closure with TPE gasket. Excellent for large batches of culture media, distilled water and other solutions.

Cat. No.	Code	Capacity, litres
BO/36310	2250-0020	10
BO/36320	2250-0050	20
BO/36350	2250-0130	50



BO/36320

Bottles, carboy, heavy-duty, PP, Autoclavable

Azlon

- Polypropylene providing good chemical resistance
- Squat form for easy location on shelves, bench top and in transit
- Moulded in handle for safe and secure carrying
- Moulded graduations for quick easy volume reference
- Autoclavable for easy cleaning and sterilisation
- Caps have moulded drill guides for easy location of probes and feed-lines.

Cat. No.	Code	Capacity, litres	Cap size, mm	Height with ca px max width mm
BO/16905	BNP05B	5	45	270 X 220
BO/16910	BNP10B	10	90	360 X 280
BO/16920	BNP20B	20	90	420 X 350



BO/169