



BB/9332

Brushes, burette

Bristle head with wire handle. Packs of 10.

Cat. No.	Head dia. mm	Overall Length, mm	For burette ml
BB/93322	13	500	10
BB/93324	16	600	25
BB/93326	19	750	50
BB/93328	25	900	100



BB/93334

Brushes, cylinder

Nylon head with fan shaped end and wire handle. Packs of 5.

Cat. No.	Head dia. X length, mm	Overall Length, mm
BB/93332	65 x 120	650
BB/93334	65 x 140	650



BB/9371

Brushes, dabber

Soft hair bristles in metal mount with wooden handle. Each.

Cat. No.	Head dia. mm	Overall Length, mm
BB/93712	9	30
BB/93714	11	34
BB/93716	14	38



BB/93340

Brushes, flask

Bristle head with hinged central joint and rubber covered stem with wooden handle. Packs of 5.

Cat. No.	Head dia. X length, mm	Overall Length, mm
BB/93340	50 x 100 3	80



BB/93752

Brushes, flat form

Hair bristles in metal mount on wooden handle. Suitable for general dusting.

Cat. No.	Brush width
BB/93752	25mm