

Burettes, PTFE key, class 'AS'

PYREX®

- Complies with ISO 385
- Manufactured from precision-bore tubing (borosilicate glass) for high accuracy
- Each burette individually engraved with serial number for complete traceability (certified burettes have a 5 point certificate)
- Each burette individually engraved with outflow time to confirm performance
- Graduations in white enamel for increased visibility with dark solutions
- Grease-free interchangeable PTFE key

Cat. No.	Capacity, ml	Grads., ml	Tol ± ml
Uncertified			
BU/12402	10	0.02	0.02
BU/12404	25	0.1	0.05
BU/12406	50	0.1	0.05
BU/12408	100	0.2	0.1
With certificate			
BU/12412	10	0.02	0.02
BU/12414	25	0.1	0.05
BU/12416	50	0.1	0.05
BU/12418	100	0.2	0.1




BU/1240

Burettes, class 'AS', borosilicate glass, P.T.F.E. key

volac

Single bore with interchangeable P.T.F.E. key and white graduations, complies with BS 846. Available uncertified or with a works certificate.

Cat. No.	Code	Capacity, ml
Uncertified		
BU/52022	RB9806/NC/Y	10 x 0.02
BU/52024	RB9806/NC/G	25 x 0.1
BU/52026	RB9806/NC/H	50 x 0.1
BU/52028	RB9806/NC/I	100 x 0.2
With certificate		
BU/52032	RB9806/WAC/Y	10 x 0.02
BU/52034	RB9806/WAC/G	25 x 0.1
BU/52036	RB9806/WAC/H	50 x 0.1
BU/52038	RB9806/WAC/I	100 x 0.2



BU/5202

Burettes, class 'B', Rotaflo®, P.T.F.E. key

MBL®

Soda glass, single bore with amber graduations. Rotaflo® stopcocks require no lubrication, do not seize with caustic solutions and give fine drop control, complies with BS846.

Cat. No.	Capacity, ml
BU/12202	10 x 0.1*
BU/12204	25 x 0.1
BU/12206	50 x 0.1
BU/12208	100 x 0.2

*supplementary to BS range.



BU/1220