



BU/5205

Burettes, double bore

Burettes, double bore, class 'B' glass key.

Soda glass with amber graduations, complies with BS 846.

Cat. No.	Capacity, ml
BU/52050	10 x 0.1
BU/52052	25 x 0.1
BU/52054	50 x 0.1
BU/52056	100 x 0.2

Burettes, automatic zero

Burettes, Sofnol/Nalfloc, water testing

Soda glass automatic zero, with reservoir

Cat. No.	Capacity, ml
BU/12006	10 x 0.1
BU/12008	20 x 0.1

Replacement/spare parts are available, please ask for details.



BU/12006

Burettes, automatic zero, Class 'AS'

PYREX®

Borosilicate glass, with PTFE key and white ceramic graduations. Burettes and reservoir bottles should be ordered separately.

Cat. No.	Capacity, ml	Joint size	Tol ± ml
Burettes			
BU/12010	10 x 0.02	29/32	0.02
BU/12012	25 x 0.05	29/32	0.03
BU/12016	50 x 0.10	29/32	0.05

Cat. No.	Code	Description
Accessories		
BU/12020	3297/12	Reservoir, glass, 2000ml, 29/32 ground neck
BU/12022	3297/14AS	Rubber blowing bulb

Burettes, Digital

Burette, digital, Digitrate Pro

Digitrate Pro™ dispenser for precise titrations with an accuracy of $\pm 0.2\%$, the head is rotatable through 360° with a see through barrel providing bubble-free titration. The borosilicate glass barrel and PTFE plunger makes it possible to use this unit with aqueous solutions and some aggressive liquids.

The scale is divided in 0.01ml divisions with an inaccuracy of $\leq 1\%$

- Continuous display of volume dispensed.
- Zero reset feature for multiple titrations
- Automatic switch-off when left unattended, preserving previous titration volume
- All wetted parts are chemically resistant and autoclavable at $+121^\circ\text{C}$
- Safety pedestal valve with off function to prevent accidental delivery
- Borosilicate glass barrel, protected with a transparent PP protection sleeve, easily removed for cleaning
- 3.6V Lithium battery giving minimum life of 60,000 three-minute titrations
- Low battery warning
- Simple user calibration and reset for ISO 9000 requirements
- Standard 30mm neck provided with 3 adapters 33, 38 and 45mm

Cat. No.	Description
BD/15006	Burette, digital 0 to 30ml x 0.01ml
BD/15008	Burette, digital 0 to 50ml x 0.01ml

Accessories

Reservoirs, amber borosilicate glass with side neck for refilling without removing the burette.

Cat. No.	Description
BD/12002	1400
BD/12004	2500



BD/15008