

Aspirators

Aspirators, borosilicate glass

PYREX®

With ground side socket (size 29/32) for stopcock assembly and unground neck to accept rubber stopper.

Cat. No.	Code	Capacity, litres	Height x dia., mm	Inlet neck dia., mm
AS/12002	1541/12D	2	248 x 138	28
AS/12004	1541/16D	5	318 x 181	44
AS/12006	1541/18D	10	298 x 227	59
AS/12008	1541/22D	20	492 x 288	59
Stopcocks				
AS/12012	1571/02M	Traditional glass key Pyrex, 29/32 4mm bore		
AS/12010	1568/04M	Rotaflo, 29/32 3mm bore		



AS/12008

Aspirators, borosilicate glass

With ground inlet and glass stopper at top and ground outlet with stopcock near base.

Cat. No.	Capacity litres	Outlet size.	Inlet neck size
AS/62002	2	19/26	24/29
AS/62005	5	19/26	29/32
AS/62010	10	24/29	45/40
AS/62020	20	29/32	50/42



AS/62005

Aspirators with stopcock, HDPE narrow mouth

Round with heavy duty carrying handles

Cat. No.	Capacity, litres
AS/56102	2.5
AS/56105	5
AS/56110	10
AS/56120	20
AS/56125	25
Spare stopcock	
AS/56101	Stopcock for 2.5, 5 and 10 litres
AS/56103	Stopcock for 20 and 25 litres



AS/561

Aspirator bottles, carboy with stopcock, LDPE

NALGENE®

With polypropylene stopcock and screw closure, size 83B.

Cat. No.	Code	Capacity, litres
AS/36304	2318-0010	4
AS/36310	2318-0020	10
AS/36320	2318-0050	20
AS/36325	2318-0065	25
AS/36350	2318-0130	50



AS/363