



BB/9471

### Burner, Meker pattern

Available in two sizes and for use with either natural (Methane) gas or LPG. The top burner is fitted with a perforated steel head which gives off multiple small but very hot flames. Each burner is supplied with 4 jets.

Cat. No.	Type	O.D. of burner head	Height
BB/94712	Natural Gas	25mm	133mm
BB/94714	L.P.G.	25mm	133mm
BB/94716	Natural Gas	43mm	180mm
BB/94718	L.P.G.	43mm	180mm



BB/94708

### Burner, Cartridge

Burner, brass with a replaceable plated steel perforated head, for use with Propane/Butane mix gas cartridges. When used with a 220gm cartridge the typical burn time is 2.5 hours. The replacement gas cartridges have been designed to enable the burner to be detached and reattached as often as needed.

Available with or without cylinders which are sold singly.

Cat. No.	Description
BB/94706	Burner head only
BB/94708	Burner head and gas cylinder complete
BB/94710	Propane/ Butane gas mix cartridge only




BB/931

### Burner, spirit, safety

Glass spirit burner, 100ml with metal cap and screw top wick holder to prevent spillage.

Cat. No.	Description
BB/93110	Spirit burner complete
BB/93112	Spare wick, 1 meter length, for cutting to size



BB/980

### Mats, ceramic, asbestos free

Cat. No.	Size, mm	Pk Qty
BB/98002	150 x 150	10
BB/98004	230 x 230	5
BB/98006	300 x 300	5
BB/98008	375 x 375	1
BB/98010	450 x 450	1



BB/93004

### Gas lighter, triggered spark

Piezo electric

Cat. No.	Description
BB/93004	Gas Lighter